### 1 Safety

The TD Pilot device has been tested and approved as compliant to all the Specifications and Standards listed in 2 *Compliance Information*, page 3 of this manual and in the 6 *Technical Specifications*, page 4 -- including, but not limited to, the Medical Device Standard (Class 1/Type B). Nevertheless, in order to ensure safe operation of your TD Pilot device, there are a few safety warnings to bear in mind:

Any serious incident that has occurred in relation to the TD Pilot should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**No modification of this equipment is allowed.**

Repairs to a Tobii Dynavox device must only be performed by Tobii Dynavox or a Tobii Dynavox authorized and approved repair center.

**Contraindication:** The TD Pilot device should never be, for the user, the only means of communicating important information. In case of failure of the TD Pilot device, the user cannot communicate using it.

The TD Pilot is not to be used as a life supporting device, and it shall not be relied on in case of loss of function due to power loss or other causes.

There could be a choking hazard risk if small parts detach from the TD Pilot device.

The TD Pilot device shall not be exposed to or used in weather conditions outside the Technical Specification of the TD Pilot device.

The TD Pilot device shall only be used with TD Pilot specific Accessories that have a mounting instruction included with the accessory.

The charging cable could present strangulation hazards to young children. Never leave small children unattended with the charging cable.

Young children or people with cognitive disabilities should not have access to, or the use of, the TD Pilot device without parental or guardian supervision.

In case the TD Pilot device malfunctions or an ESD event, restart the device.

Do not attach any decorations, stickers, papers, or similar on the screen side of the TD Pilot device. These may interfere with the Eye tracking or touch screen performance.

Never force a connector into a port. If the connector and port do not join with reasonable ease, they probably do not match. Make sure that the connector matches the port and that you have positioned the connector correctly in relation to the port.

#### 1.1 Avoiding Hearing Damage

Permanent hearing loss may occur if earphones, headphones or speakers are used at high volume. To prevent this, the volume should be set to a safe level. You can become desensitized over time to high sound levels which may then sound acceptable yet still could damage your hearing. If you experience symptoms such as ringing in your ears, please lower the volume or stop using the earphones/headphones. The louder the volume, the less time is required before your hearing could be affected.

Hearing experts suggest the following measures to protect your hearing:

- Limit the amount of time you use earphones or headphones at high volume.
- Avoid turning up the volume to block out noisy surroundings.
- Turn the volume down if you cannot hear people speaking near you.

To establish a safe volume level:

- Set your volume control at a low setting.
- Slowly increase the sound until you can hear it comfortably and clearly, without distortion.

The TD Pilot device can produce sounds in decibel ranges that may cause hearing loss for a normal hearing person, even when exposed to for less than a minute. The maximum sound level of the unit is in parity with the sound levels that a healthy young person can produce while screaming. Since the TD Pilot device is intended as a voice prosthetic, it shares the same possibilities and potential risks of causing harm to hearing. The higher decibel ranges are offered to enable communication in a noisy environment and should be used with care and only when needed in noisy environments.

## 1.2 Power Supply and Batteries

The TD Pilot device uses lithium-ion batteries. These batteries have a storage temperature range of between -20 °C/-4 °F and 40 °C/104 °F within 3 months.

Move the TD Pilot device and the batteries to a cooler environment to let the batteries charge properly.

Avoid exposing the batteries to fire or to temperatures above 50 °C/122 °F. These conditions may cause a battery to malfunction, generate heat, ignite or explode. Be aware that it is possible, in a worst case scenario, for temperatures to reach greater than those stated above in, for example, the trunk of a car on a hot day. So, storing the device, with batteries installed, in a hot car trunk could conceivably lead to a malfunction.

Do not disassemble or harm the battery. Follow the environmental laws and regulations which apply in your area when disposing of batteries.

The battery may only be replaced by the user with a TDBW1 battery pack sold by Tobii Dynavox. There is a risk of explosion if the battery is replaced with an incorrect type.

For safe operation of the TD Pilot device, use only the chargers, batteries and accessories approved by Tobii Dynavox.

Do not open (except for the battery cover), or modify, the casing of the TD Pilot device or of the power supply, since you may be exposed to potentially hazardous electrical voltage. The device contains no serviceable parts. If the TD Pilot device or its accessories are mechanically damaged, **do not to use them**.

If the battery is not charged or the TD Pilot is not connected to the power supply, the TD Pilot device will shut down.

If the Power Supply Cord becomes damaged, contact Tobii Dynavox for replacement.

Do not connect any devices with a non-medical grade power supply to any connector on the TD Pilot device. Furthermore, all configurations shall comply with the system standard IEC 60601-1. Anyone who connects additional equipment to the signal input part or signal output part is configuring a medical system and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1. The unit is for exclusive interconnection with IEC 60601-1 certified equipment in the patient environment and IEC 60601-1 certified equipment out-side of the patient environment. If in doubt, consult the technical services department or your local representative.

The appliance coupler of the power supply or separable plug is used as the Mains Disconnection Device, please do not position the TD Pilot device so that it is difficult to operate the disconnection device.

Special regulations apply to shipping Lithium-ion batteries. If dropped, crushed, punctured, thrown, abused or short-circuited, these batteries can release dangerous amounts of heat and may ignite, and are dangerous in fires.

Please reference IATA regulations when shipping lithium metal or lithium-ion batteries or cells: <http://www.iata.org/whatwedo/cargo/dgr/Pages/lithium-batteries.aspx>

Do not leave batteries uncharged for long periods to avoid deep discharge.

The TD Pilot should only be used with the NGE60-TD power supply.

## 1.3 Mounting

The TD Pilot should be mounted according to the instructions of the approved mounts used. Tobii Dynavox or its agents are not liable for damage or injuries to a person or its property due to a TD Pilot falling from a mounted configuration. The mounting of a TD Pilot is done entirely at the user's own risk.

## 1.4 Emergency

Do not rely on the device for emergency calls or banking transactions. We recommend having multiple ways to communicate in emergency situations. Banking transactions should only be carried out with a system recommended by, and approved according to the standards of your bank.

## 1.5 Infrared

The TD Pilot emits pulsed infrared (IR) light from the eye tracker. Other devices controlled by IR or susceptible to disturbance by IR light might be affected by the IR light emitted from the TD Pilot. Do not use the TD Pilot in the vicinity of such devices if their functionality is of critical importance.

## 1.6 Epilepsy

Some people with **Photosensitive Epilepsy** are susceptible to epileptic seizures or loss of consciousness when exposed to certain flashing lights or light patterns in everyday life. This may happen even if the person has no medical history of epilepsy or has never had any epileptic seizures.

A person with Photosensitive Epilepsy would also be likely to have problems with TV screens, some arcade games, and flickering fluorescent bulbs. Such people may have a seizure while watching certain images or patterns on a monitor, or even when exposed to the light sources of an eye tracker. It is estimated that about 3-5% of people with epilepsy have this type of Photosensitive Epilepsy. Many people with Photosensitive Epilepsy experience an "aura" or feel odd sensations before the seizure occurs. If you feel odd during use, move your eyes away from the eye tracker.

## 1.7 Electricity

Do not open the casing, except for the battery cover, of the TD Pilot device, since you may be exposed to potentially hazardous electrical voltage. The device contains no user serviceable parts.

## 1.8 Software

Software other than that which is pre-installed on the TD Pilot is installed at the user's own risk. External software could cause the TD Pilot to malfunction and might not be covered by the warranty.

## 1.9 Child Safety

The TD Pilot is an advanced computer system and electronic device. As such it is composed of numerous separate, assembled parts. In the hands of a child these parts have the possibility of being separated from the device, possibly constituting a choking hazard or another danger to the child.

Young children should not have access to, or the use of, the device without parental or guardian supervision.

## 1.10 Eye Tracking

Some people may experience a certain amount of fatigue (due to intentional eye focusing and hard concentration) or even a dryness of the eyes (due to less frequent blinking) when first getting used to Eye tracking. If you are experiencing fatigue or dry eyes start off slowly and limit the length of your Eye tracking sessions to your comfort level.

If needed, consult a healthcare professional regarding the use of re-moisturizing eye drops.

## 1.11 Third Party

Tobii Dynavox assumes no responsibility for any consequence resulting from use of the TD Pilot in a manner inconsistent with its intended use, including any use of the TD Pilot with third-party software and/or hardware that changes the intended use.

# 2 Compliance Information

The TD Pilot is CE marked in accordance with Regulation (EU) 2017/745 (MDR) and complies with applicable harmonized standards and the General Safety and Performance Requirements (GSPR).

## 2.1 FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation.

Modifications not expressly approved by Tobii Dynavox could void the user's authority to operate the equipment under FCC rules.

### 2.1.1 For P15B Equipment

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### 2.1.2 For Portable Devices

FCC RF Radiation Exposure Statement:

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This device was tested for typical hand held operations with the device contacted directly to the human body to the sides of the device. To maintain compliance with FCC RF exposure compliance requirements, avoid direct contact to the transmitting antenna during transmitting.

## 2.2 CE Statement

This product is CE marked as an accessory to a medical device in accordance with Regulation (EU) 2017/745 (MDR) and conforms to the applicable General Safety and Performance Requirements (GSPR).

## 2.3 Directives and Standards

The TD Pilot complies with the following directives:

- Medical Device Regulation (EU) 2017/745
- Low voltage Directive 2014/35/EU
- Electromagnetic Compatibility (EMC) Directive 2014/30/EU
- RoHS2 Directive 2011/65/EU
- WEEE Directive 2012/19/EU
- Reach Directive 2006/121/EC, 1907/2006/EC Annex 17
- ISO 14971:2019
- ISO 13485:2016

The TD Pilot device has been tested to comply with IEC/EN 60601-1 Ed 3.1, IEC/EN 62368-1, ISO 14971:2019 and other relevant standards for the intended markets.

## 3 Customer Support

For support, please contact your local representative or Support at Tobii Dynavox. In order to receive assistance as quickly as possible, make sure you have access to your TD Pilot device and, if possible, an internet connection. You should also be able to supply the serial number of the device, which you find under the Adjustable Stand on the TD Pilot Base.

For further product information and other support resources, please visit the Tobii Dynavox website [www.tobiidynavox.com](http://www.tobiidynavox.com).

## 4 Disposing of the Device

Do not dispose of the TD Pilot device in general household or office waste. Follow your local regulations for the disposal of electrical and electronic equipment.

## 5 Disposing of Batteries

Do not dispose of the batteries in general household or office waste. Follow your local regulations for the disposal of batteries.

## 6 Technical Specifications

### 6.1 Device

<b>Model</b>	<b>TD Pilot</b>
<b>Type</b>	TD Pilot
<b>Operating System</b>	Apple iPadOS
<b>CPU</b>	Apple M4 chip (9-core CPU) or newer
<b>Storage</b>	256 GB
<b>Screen Resolution</b>	2752 x 2064
<b>Screen Size</b>	13"
<b>Rear Display</b>	480 × 128 pixels
<b>Dimensions (WxHxD) TD Pilot</b>	30,4 × 25,5 × 9,0 cm 12.0 × 10.0 × 3.5 inches
<b>Weight TD Pilot</b>	2.11 kg 4.65 lbs
<b>Microphone</b>	1×Microphone
<b>Speakers</b>	2 × 10 W Closed Box Speakers

Model	TD Pilot
Connectors	<p>1 × Thunderbolt/USB 4 (iPadOS device)</p> <p>1 × USB-C</p> <p>2 × 3.5 mm Switch connector interface, (Pin out for mono plug: Sleeve = Common ground, Tip = Signal)</p> <p>1 × 3.5 mm" Headphone jack (stereo) with jack detection</p> <p>1×USB-C Power Connector</p>
Buttons	<p>1 × Top button (iPadOS device)</p> <p>1 × Volume up/down (iPadOS device)</p> <p>1 × Power On</p> <p>1 × Track Status</p>
Bluetooth®	Bluetooth 5.3
Eye Tracker (optional)	Tobii IS5TDL Module
Expected Service Life	5 years
Typical Average Battery Run Time	~10 hours
Battery Charge Time	Maximum 4 h
Deskstand	Integrated
Mounting Systems Supported	Tobii Dynavox QR adapter plate for Daessy and REHAdapt
Power Supply	15VDC, 3A, 45 W or 20VDC, 3A, 60 W AC Adapter
IP class	<p>IP54</p> <p>Pass for device only, with I/O covers in place.</p>
	<p>IP22</p> <p>Without I/O covers.</p>

## 6.2 Power Adapter

Item	Specification
Trademark	Tobii Dynavox
Manufacturer	MEAN WELL Enterprise Co., Ltd
Model Name	NGE60-TD
Rated Input	100-240Vac, 50/60Hz, 1.5-0.8A
Rated Output	5V/9V/12V/15V/20Vdc, 3A, 60W max
Output Plug	Compatible with USB-C Power Delivery up to 60W

## 6.3 Battery pack

Item	Specification	Remark
Battery Technology	Li-Ion rechargeable battery pack with gas gauge (SMBus v1.1 interface)	
Cell	6× NCR18650GA	

Item	Specification	Remark
<b>Battery Pack Capacity</b>	71,28 Wh	Initial capacity, new battery pack
<b>Nominal Voltage</b>	10.8 Vdc, 6600mAh	
<b>Charge Time</b>	Maximum 4 h	Charge from 10 to 90%
<b>Cycle Life</b>	300 cycles	Minimum 75% of initial capacity remaining
<b>Allowable Operating Temperature</b>	0 – 45 °C, 45-85%RH	Charge condition
	-20 – 60 °C, 45-85%RH	Discharge condition
<b>Storage Temperature</b>	-20 – 35 °C, 45-85%RH	1 year
	-20 – 40 °C, 45-85%RH	6 months
	-20 – 45 °C, 45-85%RH	1 month
	-20 – 50 °C, 45-85%RH	1 week
<b>Storage Time</b> <sup>1</sup>	Maximum 6 months @ charge ≥ 40%	 Do not store battery packs for a long time with less than 40% charge level.
<p>1. It is recommended that the battery shall not be stored in the device if the device is not going to be used within 6 months. If the battery is removed, the battery will not be drained as quickly as when stored in the device.</p>		

## 6.4 Eye Tracker

If installed

Technical Specifications	Tobii IS5TDL Module
<b>Working Distance</b>	45 — 95 cm 20 — 37 inches
<b>Freedom of Head Movement</b> <sup>1</sup> (Width x Height)	~20 × 20cm (7.9 × 7.9 inches) @ 50 cm from screen ~35 × 35 cm (13.8 × 13.8 inches) between 65 – 80 cm from screen
<b>Positioning</b> Distance (from screen) Track box size (Width × Height) Track box depth	45 — 95 cm (20 — 37 in) 20 × 20 — 35 × 35 cm (7.9 × 7.9 — 13.8 × 13.8 in) 50 cm (19.7 in)
<b>Gaze Data Rate</b>	33 Hz
<b>Gaze Sample Rate</b>	133 Hz
<b>Eye tracking technique</b>	Video-based pupil- and corneal- reflection eye tracking with dark and bright pupil illumination modes.
<b>Outdoor Usage Possible</b>	Yes
<b>User Calibration</b> (former tracking robustness)	>98%
<b>Detect Gaze</b> Interaction >30Hz	98% for 95% of population <sup>2</sup>
<b>Gaze Accuracy</b> Across 95% of population <sup>3</sup>	<1.58 degrees

Technical Specifications	Tobii IS5TDL Module
<b>Gaze Precision</b> Across 95% of population <sup>3</sup>	<0.2°
<b>Max Head Movement Speed</b> Eye Position Gaze Data	40 cm/s (15.7 in/s) 10 cm/s (3.9 in/s)
<b>Max Head Tilt</b>	25°
<b>Max Yaw, pitch</b>	25°
<b>Data Flow and Data Rate</b> Gaze Latency Gaze Recovery	17 ms 0 ms
<b>Mounting</b>	Built-in
<b>Power Supply</b>	Built-in
<ol style="list-style-type: none"> <li>Freedom of head movement describes the volume in front of the tracker in which the user must have at least one of the eyes. The numbers are specified parallel/orthogonal to screen surface.</li> <li>The test population excluded those who would wear correctional glasses with a +5.00 diopter or higher or suffer from an eye disease.</li> <li>The degree-of-accuracy and precision over percentage-of-population numbers result from extensive testing across representatives of the whole population. We have made use of hundreds of thousands of diagnostic images and performed tests on approximately 800 individuals with different conditions, vision, ethnicities, everyday dust, smudges or blemishes around their eyes, with eyes out of focus, etc. This has resulted in a much more robust and high-performance eye tracking experience and a much more realistic representation of true performance across the whole population, not just in a mathematically "ideal" scenario. The "ideal" degree numbers are the previous standard of measuring accuracy and precision, both previously from Tobii and presently from all eye tracking competitors. Though "ideal" numbers are useful to get a general feel about comparative quality and performance, they are not applicable to real world usage in the same way as quantitative degree-of-accuracy and precision over percentage-of-population numbers based on extensive testing across representatives of the whole population.</li> </ol>	

## 7 Guidance and Manufacturer's Declaration

Below cables information are provided for EMC reference

Cable	Max. cable length	Shielded/unshielded	Number	Cable classification
AC Power Cord	0.9 m	Unshielded	1 Set	AC Power
DC Power Cord	1.65 m	Shielded	1 Set	DC Power
Two Switch Button Cables	1.44 m	Shielded	1 Set	Signal
USB cable	0.26 m	Shielded	1 Set	Signal

### Important information regarding Electro Magnetic Compatibility (EMC)

This electrical medical equipment needs special precautions regarding EMC and put into service according to the EMC information provided in the user manual; The equipment conforms to this IEC 60601-1-2:2014+A1:2020 standard for both immunity and emissions. Nevertheless, special precautions need to be observed:

- The equipment with no ESSENTIAL PERFORMANCE is intended used in Home healthcare environment.
- WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally".
- The use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TD Pilot, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result."

- **WARNING:** If the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas, before using this equipment, it should be observed to verify that it is operating normally to assure that the equipment remains safe with regard to electromagnetic disturbances throughout the expected service life.

**STATEMENT:** For the purpose of its operation, the equipment has wireless communication function, it includes RF transmitter and receiver, 2.4 GHz, Pulse modulation.

**STATEMENT:** The equipment is designed compatible with high frequency surgical equipment; the condition includes working or standby in close proximity to high frequency surgical equipment.

Table 1: EMI Compliance Table — Emission

Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

Table 2: EMS Compliance Table — Enclosure Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10 V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3: EMS Compliance Table — Proximity fields from RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Home healthcare environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ±5kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Home healthcare environment
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Table 4: EMS Compliance Table — Input a.c. power Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Electrical fast transients/burst	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15 MHz - 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 8 0MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% U <sub>T</sub> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% U <sub>T</sub> ; 1 cycle and 70% U <sub>T</sub> ; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% U <sub>T</sub> ; 250/300 cycles

Table 5: EMS Compliance Table — Signal input/output parts Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15 MHz - 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

### 1 Safety

The TD Pilot device has been tested and approved as compliant to all the Specifications and Standards listed in 2 *Compliance Information*, page 12 of this manual and in the 6 *Technical Specifications*, page 13 -- including, but not limited to, the Medical Device Standard (Class 1/Type B). Nevertheless, in order to ensure safe operation of your TD Pilot device, there are a few safety warnings to bear in mind:

Any serious incident that has occurred in relation to the TD Pilot should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

**No modification of this equipment is allowed.**

Repairs to a Tobii Dynavox device must only be performed by Tobii Dynavox or a Tobii Dynavox authorised and approved repair centre.

**Contraindication:** The TD Pilot device should never be, for the user, the only means of communicating important information. In case of failure of the TD Pilot device, the user cannot communicate using it.

The TD Pilot is not to be used as a life supporting device, and it shall not be relied on in case of loss of function due to power loss or other causes.

There could be a choking hazard risk if small parts detach from the TD Pilot device.

The TD Pilot device shall not be exposed to or used in weather conditions outside the Technical Specification of the TD Pilot device.

The TD Pilot device shall only be used with TD Pilot specific Accessories that have a mounting instruction included with the accessory.

The charging cable could present strangulation hazards to young children. Never leave small children unattended with the charging cable.

Young children or people with cognitive disabilities should not have access to, or the use of, the TD Pilot device without parental or guardian supervision.

In case the TD Pilot device malfunctions or an ESD event, restart the device.

Do not attach any decorations, stickers, papers, or similar on the screen side of the TD Pilot device. These may interfere with the Eye tracking or touch screen performance.

Never force a connector into a port. If the connector and port do not join with reasonable ease, they probably do not match. Make sure that the connector matches the port and that you have positioned the connector correctly in relation to the port.

#### 1.1 Avoiding Hearing Damage

Permanent hearing loss may occur if earphones, headphones or speakers are used at high volume. To prevent this, the volume should be set to a safe level. You can become desensitised over time to high sound levels which may then sound acceptable yet still could damage your hearing. If you experience symptoms such as ringing in your ears, please lower the volume or stop using the earphones/headphones. The louder the volume, the less time is required before your hearing could be affected.

Hearing experts suggest the following measures to protect your hearing:

- Limit the amount of time you use earphones or headphones at high volume.
- Avoid turning up the volume to block out noisy surroundings.
- Turn the volume down if you cannot hear people speaking near you.

To establish a safe volume level:

- Set your volume control at a low setting.
- Slowly increase the sound until you can hear it comfortably and clearly, without distortion.

TheTD Pilot device can produce sounds in decibel ranges that may cause hearing loss for a normal hearing person, even when exposed to for less than a minute. The maximum sound level of the unit is in parity with the sound levels that a healthy young person can produce while screaming. Since theTD Pilot device is intended as a voice prosthetic, it shares the same possibilities and potential risks of causing harm to hearing. The higher decibel ranges are offered to enable communication in a noisy environment and should be used with care and only when needed in noisy environments.

## 1.2 Power Supply and Batteries

TheTD Pilot device uses lithium-ion batteries. These batteries have a storage temperature range of between -20 °C/-4 °F and 40 °C/104 °F within 3 months.

Move theTD Pilot device and the batteries to a cooler environment to let the batteries charge properly.

Avoid exposing the batteries to fire or to temperatures above 50 °C/122 °F. These conditions may cause a battery to malfunction, generate heat, ignite or explode. Be aware that it is possible, in a worst case scenario, for temperatures to reach greater than those stated above in, for example, the trunk of a car on a hot day. So, storing the device, with batteries installed, in a hot car trunk could conceivably lead to a malfunction.

Do not disassemble or harm the battery. Follow the environmental laws and regulations which apply in your area when disposing of batteries.

The battery may only be replaced by the user with a TDBW1 battery pack sold by Tobii Dynavox. There is a risk of explosion if the battery is replaced with an incorrect type.

For safe operation of theTD Pilot device, use only the chargers, batteries and accessories approved by Tobii Dynavox.

Do not open (except for the battery cover), or modify, the casing of theTD Pilot device or of the power supply, since you may be exposed to potentially hazardous electrical voltage. The device contains no serviceable parts. If theTD Pilot device or its accessories are mechanically damaged, **do not to use them**.

If the battery is not charged or theTD Pilot is not connected to the power supply, theTD Pilot device will shut down.

If the Power Supply Cord becomes damaged, contactTobii Dynavox for replacement.

Do not connect any devices with a non-medical grade power supply to any connector on theTD Pilot device. Furthermore, all configurations shall comply with the system standard IEC 60601-1. Anyone who connects additional equipment to the signal input part or signal output part is configuring a medical system and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1. The unit is for exclusive interconnection with IEC 60601-1 certified equipment in the patient environment and IEC 60601-1 certified equipment out-side of the patient environment. If in doubt, consult the technical services department or your local representative.

The appliance coupler of the power supply or separable plug is used as the Mains Disconnection Device, please do not position theTD Pilot device so that it is difficult to operate the disconnection device.

Special regulations apply to shipping Lithium-ion batteries. If dropped, crushed, punctured, thrown, abused or short-circuited, these batteries can release dangerous amounts of heat and may ignite, and are dangerous in fires.

Please reference IATA regulations when shipping lithium metal or lithium-ion batteries or cells: <http://www.iata.org/whatwedo/cargo/dgr/Pages/lithium-batteries.aspx>

Do not leave batteries uncharged for long periods to avoid deep discharge.

TheTD Pilot should only be used with theNGE60-TD power supply.

## 1.3 Mounting

TheTD Pilot should be mounted according to the instructions of the approved mounts used.Tobii Dynavox or its agents are not liable for damage or injuries to a person or its property due to aTD Pilot falling from a mounted configuration. The mounting of aTD Pilot is done entirely at the user's own risk.

## 1.4 Emergency

Do not rely on the device for emergency calls or banking transactions. We recommend having multiple ways to communicate in emergency situations. Banking transactions should only be carried out with a system recommended by, and approved according to the standards of your bank.

## 1.5 Infrared

TheTD Pilot emits pulsed infrared (IR) light from the eye tracker. Other devices controlled by IR or susceptible to disturbance by IR light might be affected by the IR light emitted from the TD Pilot. Do not use theTD Pilot in the vicinity of such devices if their functionality is of critical importance.

## 1.6 Epilepsy

Some people with **Photosensitive Epilepsy** are susceptible to epileptic seizures or loss of consciousness when exposed to certain flashing lights or light patterns in everyday life. This may happen even if the person has no medical history of epilepsy or has never had any epileptic seizures.

A person with Photosensitive Epilepsy would also be likely to have problems with TV screens, some arcade games, and flickering fluorescent bulbs. Such people may have a seizure while watching certain images or patterns on a monitor, or even when exposed to the light sources of an eye tracker. It is estimated that about 3-5% of people with epilepsy have this type of Photosensitive Epilepsy. Many people with Photosensitive Epilepsy experience an "aura" or feel odd sensations before the seizure occurs. If you feel odd during use, move your eyes away from the eye tracker.

## 1.7 Electricity

Do not open the casing, except for the battery cover, of theTD Pilot device, since you may be exposed to potentially hazardous electrical voltage. The device contains no user serviceable parts.

## 1.8 Software

Software other than that which is pre-installed on theTD Pilot is installed at the user's own risk. External software could cause theTD Pilot to malfunction and might not be covered by the warranty.

## 1.9 Child Safety

TheTD Pilot is an advanced computer system and electronic device. As such it is composed of numerous separate, assembled parts. In the hands of a child these parts have the possibility of being separated from the device, possibly constituting a choking hazard or another danger to the child.

Young children should not have access to, or the use of, the device without parental or guardian supervision.

## 1.10 Eye Tracking

Some people may experience a certain amount of fatigue (due to intentional eye focusing and hard concentration) or even a dryness of the eyes (due to less frequent blinking) when first getting used to Eye tracking. If you are experiencing fatigue or dry eyes start off slowly and limit the length of your Eye tracking sessions to your comfort level.

If needed, consult a healthcare professional regarding the use of re-moisturising eye drops.

## 1.11 Third Party

Tobii Dynavox assumes no responsibility for any consequence resulting from use of theTD Pilot in a manner inconsistent with its intended use, including any use of theTD Pilot with third-party software and/or hardware that changes the intended use.

# 2 Compliance Information

TheTD Pilot is CE marked in accordance with Regulation (EU) 2017/745 (MDR) and complies with applicable harmonised standards and the General Safety and Performance Requirements (GSPR).

## 2.1 FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation.

Modifications not expressly approved by Tobii Dynavox could void the user's authority to operate the equipment under FCC rules.

### 2.1.1 For P15B Equipment

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### 2.1.2 For Portable Devices

FCC RF Radiation Exposure Statement:

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This device was tested for typical hand held operations with the device contacted directly to the human body to the sides of the device. To maintain compliance with FCC RF exposure compliance requirements, avoid direct contact to the transmitting antenna during transmitting.

## 2.2 CE Statement

This product is CE marked as an accessory to a medical device in accordance with Regulation (EU) 2017/745 (MDR) and conforms to the applicable General Safety and Performance Requirements (GSPR).

## 2.3 Directives and Standards

The TD Pilot complies with the following directives:

- Medical Device Regulation (EU) 2017/745
- Low voltage Directive 2014/35/EU
- Electromagnetic Compatibility (EMC) Directive 2014/30/EU
- RoHS2 Directive 2011/65/EU
- WEEE Directive 2012/19/EU
- Reach Directive 2006/121/EC, 1907/2006/EC Annex 17
- ISO 14971:2019
- ISO 13485:2016

The TD Pilot device has been tested to comply with IEC/EN 60601-1 Ed 3.1, IEC/EN 62368-1, ISO 14971:2019 and other relevant standards for the intended markets.

## 3 Customer Support

For support, please contact your local representative or Support at Tobii Dynavox. In order to receive assistance as quickly as possible, make sure you have access to your TD Pilot device and, if possible, an internet connection. You should also be able to supply the serial number of the device, which you find under the Adjustable Stand on the TD Pilot Base.

For further product information and other support resources, please visit the Tobii Dynavox website [www.tobiidynavox.com](http://www.tobiidynavox.com).

## 4 Disposing of the Device

Do not dispose of the TD Pilot device in general household or office waste. Follow your local regulations for the disposal of electrical and electronic equipment.

## 5 Disposing of Batteries

Do not dispose of the batteries in general household or office waste. Follow your local regulations for the disposal of batteries.

## 6 Technical Specifications

### 6.1 Device

<b>Model</b>	<b>TD Pilot</b>
<b>Type</b>	TD Pilot
<b>Operating System</b>	Apple iPadOS
<b>CPU</b>	Apple M4 chip (9-core CPU) or newer
<b>Storage</b>	256 GB
<b>Screen Resolution</b>	2752 x 2064
<b>Screen Size</b>	13"
<b>Rear Display</b>	480 × 128 pixels
<b>Dimensions (WxHxD) TD Pilot</b>	30.4 × 25.5 × 9.0 cm 12.0 × 10.0 × 3.5 inches
<b>Weight TD Pilot</b>	2.11 kg 4.65 lbs
<b>Microphone</b>	1×Microphone
<b>Speakers</b>	2 × 10 W Closed Box Speakers

Model	TD Pilot
Connectors	<p>1 × Thunderbolt/USB 4 (iPadOS device)</p> <p>1 × USB-C</p> <p>2 × 3.5 mm Switch connector interface, (Pin out for mono plug: Sleeve = Common ground, Tip = Signal)</p> <p>1 × 3.5 mm" Headphone jack (stereo) with jack detection</p> <p>1×USB-C Power Connector</p>
Buttons	<p>1 × Top button (iPadOS device)</p> <p>1 × Volume up/down (iPadOS device)</p> <p>1 × Power On</p> <p>1 × Track Status</p>
Bluetooth®	Bluetooth 5.3
Eye Tracker (optional)	Tobii IS5TDL Module
Expected Service Life	5 years
Typical Average Battery Run Time	~10 hours
Battery Charge Time	Maximum 4 h
Deskstand	Integrated
Mounting Systems Supported	Tobii Dynavox QR adapter plate for Daessy and REHAdapt
Power Supply	15VDC, 3A, 45 W or 20VDC, 3A, 60 W AC Adapter
IP class	<p>IP54</p> <p>Pass for device only, with I/O covers in place.</p>
	<p>IP22</p> <p>Without I/O covers.</p>

## 6.2 Power Adapter

Item	Specification
Trademark	Tobii Dynavox
Manufacturer	MEAN WELL Enterprise Co., Ltd
Model Name	NGE60-TD
Rated Input	100-240Vac, 50/60Hz, 1.5-0.8A
Rated Output	5V/9V/12V/15V/20Vdc, 3A, 60W max
Output Plug	Compatible with USB-C Power Delivery up to 60W

## 6.3 Battery pack

Item	Specification	Remark
Battery Technology	Li-Ion rechargeable battery pack with gas gauge (SMBus v1.1 interface)	
Cell	6× NCR18650GA	

Item	Specification	Remark
<b>Battery Pack Capacity</b>	71,28 Wh	Initial capacity, new battery pack
<b>Nominal Voltage</b>	10.8 Vdc, 6600mAh	
<b>Charge Time</b>	Maximum 4 h	Charge from 10 to 90%
<b>Cycle Life</b>	300 cycles	Minimum 75% of initial capacity remaining
<b>Allowable Operating Temperature</b>	0 – 45 °C, 45-85%RH	Charge condition
	-20 – 60 °C, 45-85%RH	Discharge condition
<b>Storage Temperature</b>	-20 – 35 °C, 45-85%RH	1 year
	-20 – 40 °C, 45-85%RH	6 months
	-20 – 45 °C, 45-85%RH	1 month
	-20 – 50 °C, 45-85%RH	1 week
<b>Storage Time</b> <sup>1</sup>	Maximum 6 months @ charge ≥ 40%	 Do not store battery packs for a long time with less than 40% charge level.
<p>1. It is recommended that the battery shall not be stored in the device if the device is not going to be used within 6 months. If the battery is removed, the battery will not be drained as quickly as when stored in the device.</p>		

## 6.4 Eye Tracker

If installed

Technical Specifications	Tobii IS5TDL Module
<b>Working Distance</b>	45 — 95 cm 20 — 37 inches
<b>Freedom of Head Movement</b> <sup>1</sup> (Width x Height)	~20 × 20cm (7.9 × 7.9 inches) @ 50 cm from screen ~35 × 35 cm (13.8 × 13.8 inches) between 65 – 80 cm from screen
<b>Positioning</b> Distance (from screen) Track box size (Width × Height) Track box depth	45 — 95 cm (20 — 37 in) 20 × 20 — 35 × 35 cm (7.9 × 7.9 — 13.8 × 13.8 in) 50 cm (19.7 in)
<b>Gaze Data Rate</b>	33 Hz
<b>Gaze Sample Rate</b>	133 Hz
<b>Eye tracking technique</b>	Video-based pupil- and corneal- reflection eye tracking with dark and bright pupil illumination modes.
<b>Outdoor Usage Possible</b>	Yes
<b>User Calibration</b> (former tracking robustness)	>98%
<b>Detect Gaze</b> Interaction >30Hz	98% for 95% of population <sup>2</sup>
<b>Gaze Accuracy</b> Across 95% of population <sup>3</sup>	<1.58 degrees

Technical Specifications	Tobii IS5TDL Module
<b>Gaze Precision</b> Across 95% of population <sup>3</sup>	<0.2°
<b>Max Head Movement Speed</b> Eye Position Gaze Data	40 cm/s (15.7 in/s) 10 cm/s (3.9 in/s)
<b>Max Head Tilt</b>	25°
<b>Max Yaw, pitch</b>	25°
<b>Data Flow and Data Rate</b> Gaze Latency Gaze Recovery	17 ms 0 ms
<b>Mounting</b>	Built-in
<b>Power Supply</b>	Built-in
<ol style="list-style-type: none"> <li>Freedom of head movement describes the volume in front of the tracker in which the user must have at least one of the eyes. The numbers are specified parallel/orthogonal to screen surface.</li> <li>The test population excluded those who would wear correctional glasses with a +5.00 diopter or higher or suffer from an eye disease.</li> <li>The degree-of-accuracy and precision over percentage-of-population numbers result from extensive testing across representatives of the whole population. We have made use of hundreds of thousands of diagnostic images and performed tests on approximately 800 individuals with different conditions, vision, ethnicities, everyday dust, smudges or blemishes around their eyes, with eyes out of focus, etc. This has resulted in a much more robust and high-performance eye tracking experience and a much more realistic representation of true performance across the whole population, not just in a mathematically "ideal" scenario. The "ideal" degree numbers are the previous standard of measuring accuracy and precision, both previously from Tobii and presently from all eye tracking competitors. Though "ideal" numbers are useful to get a general feel about comparative quality and performance, they are not applicable to real world usage in the same way as quantitative degree-of-accuracy and precision over percentage-of-population numbers based on extensive testing across representatives of the whole population.</li> </ol>	

## 7 Guidance and Manufacturer's Declaration

Below cables information are provided for EMC reference

Cable	Max. cable length	Shielded/unshielded	Number	Cable classification
AC Power Cord	0.9 m	Unshielded	1 Set	AC Power
DC Power Cord	1.65 m	Shielded	1 Set	DC Power
Two Switch Button Cables	1.44 m	Shielded	1 Set	Signal
USB cable	0.26 m	Shielded	1 Set	Signal

### Important information regarding Electro Magnetic Compatibility (EMC)

This electrical medical equipment needs special precautions regarding EMC and put into service according to the EMC information provided in the user manual; The equipment conforms to this IEC 60601-1-2:2014+A1:2020 standard for both immunity and emissions. Nevertheless, special precautions need to be observed:

- The equipment with no ESSENTIAL PERFORMANCE is intended used in Home healthcare environment.
- WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally".
- The use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TD Pilot, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result."

- **WARNING:** If the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas, before using this equipment, it should be observed to verify that it is operating normally to assure that the equipment remains safe with regard to electromagnetic disturbances throughout the expected service life.

**STATEMENT:** For the purpose of its operation, the equipment has wireless communication function, it includes RF transmitter and receiver, 2.4 GHz, Pulse modulation.

**STATEMENT:** The equipment is designed compatible with high frequency surgical equipment; the condition includes working or standby in close proximity to high frequency surgical equipment.

Table 1: EMI Compliance Table — Emission

Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

Table 2: EMS Compliance Table — Enclosure Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10 V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3: EMS Compliance Table — Proximity fields from RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Home healthcare environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ±5kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Home healthcare environment
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Table 4: EMS Compliance Table — Input a.c. power Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Electrical fast transients/burst	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15 MHz - 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 8 0MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% U <sub>T</sub> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% U <sub>T</sub> ; 1 cycle and 70% U <sub>T</sub> ; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% U <sub>T</sub> ; 250/300 cycles

Table 5: EMS Compliance Table — Signal input/output parts Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home healthcare environment
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15 MHz - 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

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